

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property  
Organization  
International Bureau



(43) International Publication Date  
27 May 2004 (27.05.2004)

PCT

(10) International Publication Number  
**WO 2004/044247 A2**

- (51) International Patent Classification<sup>7</sup>: **C12Q 1/70**
- (21) International Application Number: PCT/IB2003/005312
- (22) International Filing Date: 3 November 2003 (03.11.2003)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data: 60/425,327 12 November 2002 (12.11.2002) US
- (71) Applicant (*for all designated States except US*): **GENO-LIFE** [FR/FR]; Biopole Clermont Limagne, Saint Beaulieu, F-63360 Gerzat (FR).
- (72) Inventors; and
- (75) Inventors/Applicants (*for US only*): **CHAUBRON, Franck** [FR/FR]; 6, Allée de la Fontoune, F-63400 Chamalières (FR). **MARTIN-MINVIELLE, Anne, Céline** [FR/FR]; 61, rue des Gravouses, F-63400 Chamalières (FR). **GROULON, Sophie** [FR/FR]; 1, Impasse Edouard Branly, F-63800 Courmon D'Auvergne (FR).
- (74) Agents: **MARTIN, Jean-Jacques** et al.; Cabinet Regimbeau, 20, rue de Chazelles, F-75847 Paris Cedex 17 (FR).
- (81) Designated States (*national*): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.
- (84) Designated States (*regional*): ARIPO patent (BW, GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).
- Published:**  
— *without international search report and to be republished upon receipt of that report*
- For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.*

(54) Title: ONE STEP REAL-TIME RT PCR KITS FOR THE UNIVERSAL DETECTION OF ORGANISMS IN INDUSTRIAL PRODUCTS

(57) Abstract: This invention is related to a novel sample preparation, probes, couple primer sets for one step real time reverse transcriptase polymerise chain reaction (RT-PCR), methods and kits for the universal detection of alive bacteria and/or fungus-yeast in pharmaceutical, cosmetic and non clinical samples.



WO 2004/044247 A2

**One step real-time RT PCR kits for the universal detection of organisms in  
industrial products**

**TECHNICAL FIELD**

5

The invention pertains to the field of methods and reagents for detecting bacteria and fungus-yeast found in pharmaceutical, cosmetic and non clinical samples.

More specifically, the present invention relates to a sample preparation, primer sets, probe sets and methods for one step real-time RT PCR kits for the universal detection  
10 of alive bacteria and fungus-yeast in sterile or non sterile industrial product in less than 24 hours.

**BACKGROUND OF THE INVENTION**

15 The use of specific polynucleotide sequences or Peptide Nucleic Acid as primers and/or probes for the recognition of contaminant and infectious agents is becoming a valuable alternative to problematic growth requirements assays, visible (colony) growth features, microscopic morphology, staining reactions, and biochemical characteristics. Most of the time these technologies are too slow for a real used in  
20 industrial controls.

For example, PCT publication W084/02721, published Jul. 19, 1984 describes the use of nucleic acid probes complementary to targeted nucleic acid sequences composed of ribosomal RNA, transfer RNA, or other RNA in hybridization  
25 procedures to detect the target nucleic acid sequence. While the assay may provide greater sensitivity and specificity than known DNA hybridization assays, hybridization procedures which require the use of a complementary probe are

generally dependent upon the cultivation of a test organism and are, therefore, unsuitable for rapid analysis.

These probes are useful in hybridizing to RNA amplified by the Reverse Transcriptase Polymerase Chain Reaction (RT PCR). RT-PCR is a powerful  
5 ribonucleic acid amplification technique that can be used for the detection of small numbers of ribonucleotide acid targets from bacteria and/or from fungus-yeast whose in vitro cultivation is difficult or lengthy. RT PCR requires the presence of living specimens for detection. In its simplest form, RT PCR is an in vitro method for the enzymatic synthesis of specific cDNA sequences. Using one oligonucleotide primers  
10 that hybridize to RNA strand and flank the region of interest in the target cDNA, several cDNA are synthesised by Reverse Transcriptase. A repetitive series of cycles involving template denaturation, primer annealing, and extension of annealed primers by DNA polymerase results in the exponential accumulation of a specific fragment whose termini are defined by the 5' ends of the primers. PCR produce a selective  
15 enrichment of a specific DNA sequence by a factor of 10<sup>sup.12</sup>. The PCR method is described in Saiki et al, 1985, Science 230:1350 and is the subject of U.S. Pat. Nos. 4,683,195, 4,683,202 and 4,800,159 (these references are incorporated herein by reference). This method has been used to detect the presence of the aberrant sequence in the beta-globin gene which is related to sickle cell anemia (Saiki et al., 1985,  
20 supra) and the human immunodeficiency virus (HIV) RNA (Byrne et al., 1988, Nuc. Acids Res. 16:4165).

In order to successfully treat a contamination caused by a bacteria or fungus-yeast in a sterile or non sterile product from industry, a rapid and accurate detection is  
25 required. Bacterial and fungus-yeast detection have traditionally been accomplished by pure culture isolation, followed by identification procedures that make use of knowledge of specimen source, growth requirements, visible (colony) growth

features, microscopic morphology, staining reactions, and biochemical characteristics.

It is clear that a rapid diagnostic method, less than 24 hours, for detecting bacteria and fungus-yeast in industrial samples with the same sensitivity as culture would be a  
5 significant improvement over currently used methods.

#### SUMMARY OF THE INVENTION

The present invention pertains to methods and reagents for the rapid detection of  
10 bacteria and fungus --yeast in sterile and non sterile product in less than 24 hours.

In a preferred embodiment, a target region from a one-step Reverse Transcriptase Polymerase Chain Reaction of RNA and the resultant amplified DNA is treated with probes which can hybridize to the amplified DNA of bacteria or fungus-yeast but not other organisms (mammalian, plant, insects...) or virus .

15 The Tth DNA polymerase is a thermostable enzyme with RNA-dependent Reverse Transcriptase activity and with DNA-dependent Polymerase activity, allowing the combination of RT and PCR in a single-tube reaction resulting in a faster analysis of presence of RNA from bacteria, fungus-yeast.

Using one-step Real-time Reverse Transcriptase Polymerase Chain Reaction, the  
20 invention enable the user to perform a rapid RT-PCR and simultaneously detect and quantify the presence of RNA from bacteria and/or fungus-yeast by monitoring fluorescence during real time polymerase chain reaction amplification with any risk of false positive due to opening tube between RT and PCR and from possible PCR product environmental contamination due to precedent amplification reactions in the  
25 laboratory.

## DETAILED DESCRIPTION OF THE INVENTION

The methods of the present invention thus enable determination of the presence of bacteria and/or fungus-yeast more rapidly than technologies with prior art detection  
5 methods.

Using one-step Real-time Reverse Transcriptase Polymerase Chain Reaction, the invention enable the user to perform a rapid RT-PCR and simultaneously analyse and quantify the presence of RNA from bacteria and/or fungus-yeast by monitoring fluorescence during real time polymerase chain reaction amplification with any risk  
10 of false positive due to opening tube between RT and PCR and from possible PCR product environmental contamination due to precedent amplification reactions.

The basic RT PCR process is carried out as follows.

A sample is provided which needs to be tested or which is suspected of contain a  
15 particular ribonucleic acid sequence of interest, the "target sequence." The ribonucleic acid contained in the sample may be first reverse transcribed into cDNA (using enzyme like Tth DNA polymerase as purified enzyme and a oligonucleotide or PNA ), and then denatured, using physical means, which are known to those of skill in the art. A preferred physical means for strand separation involves heating the  
20 nucleic acid until it is completely (>99%) denatured. Methods for the amplification of RNA targets using a thermostable DNA polymerase are described in PCT/US90/07641, filed Dec. 21, 1990, and incorporated herein by reference.

The denatured DNA strands are then incubated in the same tube with the selected  
25 oligonucleotide primers under hybridization conditions, conditions which enable the binding of the primers to the single DNA strands. As known in the art, the primers are selected so that their relative positions along a duplex sequence are such that an extension product synthesized from one primer, when it is separated from its

complement, serves as a template for the extension of the other primer to yield a replicate chain of defined length.

The primer must be sufficiently long to prime the synthesis of extension products in the presence of the agent for polymerization. The exact length of the primers will  
 5 depend on many factors, including temperature, source of the primer and use of the method.

Preferred oligonucleotide primers for use in the present invention are selected from the group consisting of

10	Seq ID No 1	TGGAGCATGTGGTTTAATTCGA	[primer forward]
	Seq ID No 2	TGCGGGACTTAACCCAACA	[primer reverse]
	Seq ID No 3	AGAGTTTGATCATGGCTCAGA	[primer forward]
	Seq ID No 4	TTACCCACCTACTAGCTAAT	[primer reverse]
	Seq ID No 5	GYGGAGCATGTGGYTTAATTCG	[primer forward]
15	Seq ID No 6	TTGCGCTCGTTTCGGGACTT	[primer reverse]
	Seq ID No 7	GGGAAACTCACCAGGTCCA	[primer forward]
	Seq ID No 8	CGTTATCGCAATTAAGCAGACA	[primer reverse]
	Seq ID No 9	GGTAACGGGGAATWAGGGTTC	[primer forward]
	Seq ID No 10	TTGGGTAATTTGCGCGCCTG	[primer reverse]

20

Universal code : Y= (C/T), R = (A/G), W = (A/T)

Template-dependent extension of the oligonucleotide primer(s) is then catalyzed by the polymerizing agent ( in the presence of adequate amounts of the four  
 25 deoxyribonucleoside triphosphates (dATP, dGTP, dCTP, and dTTP) or analogs), in a reaction medium which is comprised of the appropriate salts, metal cations, and pH buffering system.

The products of the synthesis are duplex molecules consisting of the template strands and the primer extension strands, which include the target sequence. These products, in turn, serve as templates for another round of replication. In the second round of replication, the primer extension strand of the first cycle is annealed with its  
5 complementary primer; synthesis yields a "short" product which is bounded on both the 5'-and the 3'-ends by primer sequences or their complements. Repeated cycles of denaturation, primer annealing, and extension result in the exponential accumulation of the target region defined by the primers. Sufficient cycles are run to achieve the desired amount of polynucleotide containing the target region of nucleic acid. The  
10 desired amount may vary, and is determined by the function which the product polynucleotide is to serve.

The PCR method is performed in a fashion where all of the reagents are added simultaneously, in one step.  
15 In a preferred method, the RT PCR reaction is carried out as an automated process which utilizes a thermostable enzyme like Tth.

The types of machines used are commercially available from Roche Diagnostics (LigthCycler), Cepheid (Smart Cycler, GeneXpert), BioRad (Icycler), Corbett  
20 Research (Rotor-Gene)... and most suitable equipment developed for real time PCR assays and commercial use.

Those skilled in the art will also be aware of the problems of contamination of a PCR by the nucleic acid from bacteria previously present in water used for buffer and  
25 resulting in non specific amplification. Methods to reduce these problems are provided by using adequate filtration systems to avoid DNA strand fragments with a size higher than 100 bp. All reagents used in the RT PCR reaction have to be processed before using.

During amplification by PCR, the target polynucleotides may be detected directly by hybridization with a probe polynucleotide which forms a stable hybrid with the target sequence under high stringency to low stringency hybridization and washing conditions.

5

Probes are typically labeled with non-radioactive labeling systems, such as fluoresceins and derivated systems.

Therefore, in one embodiment, the invention relates to a method and kit for  
10 determining the presence of bacteria or fungus-yeast ribonucleic acid (RNA) in a sample suspected of containing said bacteria and/or fungus, wherein said polynucleotide comprises a selected target region, said method comprising:

(a) extract bacteria or fungus-yeast ribonucleic acid (RNA) from the sample up to  
15 1000 ml by centrifiltration on membranes and/or DEAE resin following by incubation with DNase.

(b) incubating the bacteria or fungus-yeast ribonucleic acid (RNA) with a  
thermostable enzyme with RNA-dependent Reverse Transcriptase activity and with  
20 DNA-dependent Polymerase activity, allowing the combination of RT and PCR in a single-tube reaction, such as Tth DNA polymerase or an enzyme like Tth DNA polymerase, and polynucleotide primers with a nucleotide sequence selected from the group consisting of

Seq ID No 2	TGCGGGACTTAACCCAACA	[primer reverse]
25 Seq ID No 4	TTACCCACCTACTAGCTAAT	[primer reverse]
Seq ID No 6	TTGCGCTCGTTTRCGGGACTT	[primer reverse]
Seq ID No 8	CGTTATCGCAATTAAGCAGACA	[primer reverse]
Seq ID No 10	TTGGGTAATTTGCGCGCCTG	[primer reverse]



under conditions which allow hybridization of the polynucleotide to the ribonucleotide target region and Reverse Transcriptase activity of the said polymerase, or enzyme like Tth, for cDNA synthesis; and

5

(c) amplified the cDNAs formed to a detectable level by Polymerase Chain Reaction with said polymerase enzyme like Tth DNA polymerase and polynucleotide primers and probes with a nucleotide sequence selected from the group consisting of

	Seq ID No 1	TGGAGCATGTGGTTTAATTCGA	[primer forward]
10	Seq ID No 2	TGCGGGACTTAACCCAACA	[primer reverse]
	Seq ID No 3	AGAGTTTGATCATGGCTCAGA	[primer forward]
	Seq ID No 4	TTACCCACCTACTAGCTAAT	[primer reverse]
	Seq ID No 5	GYGGAGCATGTGGYTAAATTCG	[primer forward]
	Seq ID No 6	TTGCGCTCGTTTRCGGGACTT	[primer reverse]
15	Seq ID No 7	GGGAAACTCACCAGGTCCA	[primer forward]
	Seq ID No 8	CGTTATCGCAATTAAGCAGACA	[primer reverse]
	Seq ID No 9	GGTAACGGGGAATWAGGGTTC	[primer forward]
	Seq ID No 10	TTGGGTAATTTGCGCGCCTG	[primer reverse]
	Seq ID No 11	TGCATGGYTGTCTGTCAGCTCGTG	[probe forward]
20	Seq ID No 12	GAGTGGCGGACGGGTGAGTAA	[probe forward]
	Seq ID No 13	ACAGGTGGTGCATGGTTGTC	[probe forward]
	Seq ID No 14	TCAGCTCGTGTCGTGAGATGTT	[probe forward]
	Seq ID No 15	ACAGGTGCTGCATGGCTGTC	[probe forward]
	Seq ID No 16	TCAGCTCGTGTTGTGAAATGTT	[probe forward]
25	Seq ID No 17	AGGATTGACAGATTGAGAGCTCTT	[probe forward]
	Seq ID No 18	CGGAGAGGGAGCCTGAGAA	[probe forward]
	Seq ID No 19	CGGCTACCACATCCAAGGAA	[probe forward]

The cDNA target sequence can be synthesised by Reverse Transcriptase activity of Tth or an enzyme like Tth and is amplified by the DNA-dependent Polymerase activity of the DNA polymerase in the same tube by means of one step real time RT - PCR.

5

More particularly, the composition for detecting bacteria comprises a polynucleotide primers and a probe consisting of the sequence

Seq ID No 1	TGGAGCATGTGGTTTAATTCGA	[primer forward]
Seq ID No 2	TGCGGGACTTAACCCAACA	[primer reverse]
10 Seq ID No 11	TGCATGGYTGTCGTCAGCTCGTG	[probe forward]

Alternatively, the composition for detecting bacteria comprises a polynucleotide primers and a probe consisting of the sequence

Seq ID No 3	AGAGTTTGATCATGGCTCAGA	[primer forward]
15 Seq ID No 4	TTACCCACCTACTAGCTAAT	[primer reverse]
Seq ID No 12	GAGTGGCGGACGGGTGAGTAA	[probe forward]

The invention may also be practiced with a composition for detecting bacteria which comprises a polynucleotide primers and a probe consisting of the sequence

20 Seq ID No 5	GYGGAGCATGTGGYTAAATTCG	[primer forward]
Seq ID No 6	TTGCGCTCGTTTCGGGACTT	[primer reverse]
Seq ID No 13	ACAGGTGGTGCATGGTTGTC	[probe forward]
Seq ID No 14	TCAGCTCGTGTCGTGAGATGTT	[probe forward]
Seq ID No 15	ACAGGTGCTGCATGGCTGTC	[probe forward]
25 Seq ID No 16	TCAGCTCGTGTTGTGAAATGTT	[probe forward]

The invention also concerns the above mentioned method and kit wherein the composition for detecting fungus-yeast comprises a polynucleotide primers and a probe consisting of the sequence

Seq ID No 7 GGGAAACTCACCAGGTCCA [primer forward]

5 Seq ID No 8 CGTTATCGCAATTAAGCAGACA [primer reverse]

Seq ID No 17 AGGATTGACAGATTGAGAGCTCTT [probe forward]

Alternatively, the composition for detecting fungus-yeast comprises a polynucleotide primers and a probe consisting of the sequence

10 Seq ID No 9 GGTAACGGGGAATWAGGGTTC [primer forward]

Seq ID No 10 TTGGGTAATTTGCGCGCCTG [primer reverse]

Seq ID No 18 CGGAGAGGGAGCCTGAGAA [probe forward]

Seq ID No 19 CGGCTACCACATCCAAGGAA [probe forward]

15 The preferred combination of primers and probes used for detection all bacteria and/or fungus-yeast consists of the sequence :

Seq ID No 1+ Seq ID No 2 +Seq ID No 11

or

Seq ID No 3+ Seq ID No 4 +Seq ID No 12

20 or

Seq ID No 5+ Seq ID No 6 +Seq ID No 13 + Seq ID No 14 + Seq ID No 15 +Seq ID No 16

or

Seq ID No 7+ Seq ID No 8 +Seq ID No 17

25 or

Seq ID No 9+ Seq ID No 10 +Seq ID No 18 + Seq ID No 19

or

Seq ID No 1+ Seq ID No 2 +Seq ID No 11 + Seq ID No 7+ Seq ID No 8 +Seq ID No 17

or

Seq ID No 3+ Seq ID No 4 +Seq ID No 12 + Seq ID No 7+ Seq ID No 8 +Seq ID No 17

or

Seq ID No 5+ Seq ID No 6 +Seq ID No 13 + Seq ID No 14 + Seq ID No 15 +Seq ID No 16 + Seq ID No 9+ Seq ID No 10 +Seq ID No 18 + Seq ID No 19

- 10 As mentioned above, the polynucleotide primers and probes may be natural nucleic acid or Peptide Nucleic Acid (PNA) which can hybridize to nucleic acid (DNA and RNA).

The RNA may also be quantified and compared with quantified external standard RNA from by exemple *Escherichia coli* and *Candida spp.*

15

By way of further specificity, the following probe nucleotide base pair data is provided. Preferred oligonucleotide probes for use in the present invention are selected from the group consisting of

- |    |              |                          |                 |
|----|--------------|--------------------------|-----------------|
|    | Seq ID No 11 | TGCATGGYTGTCGTCAGCTCGTG  | [probe forward] |
| 20 | Seq ID No 12 | GAGTGGCGGACGGGTGAGTAA    | [probe forward] |
|    | Seq ID No 13 | ACAGGTGGTGCATGGTTGTC     | [probe forward] |
|    | Seq ID No 14 | TCAGCTCGTGTCTGTGAGATGTT  | [probe forward] |
|    | Seq ID No 15 | ACAGGTGCTGCATGGCTGTC     | [probe forward] |
|    | Seq ID No 16 | TCAGCTCGTGTGTTGTGAAATGTT | [probe forward] |
| 25 | Seq ID No 17 | AGGATTGACAGATTGAGAGCTCTT | [probe forward] |
|    | Seq ID No 18 | CGGAGAGGGAGCCTGAGAA      | [probe forward] |
|    | Seq ID No 19 | CGGCTACCACATCCAAGGAA     | [probe forward] |

Reverse probes are not usable because in a one-step RT-PCR the probe should not hybridize the RNA sequence zone where the cDNA is synthesised by Reverse Transcriptase.

- 5 The sequence of the preferred oligonucleotide primers and probes of the invention are based on the rRNA gene. Oligonucleotide rRNA gene for the detection of nucleic acids from various microorganisms have been described in the scientific literature. For example, universal bacterial probes have been described by Wilson et al., 1990, J. Clinical Microbiology 28:1942-1946, and Chem et al., 1989, FEMS Microbiology Letters 57:19-24.
- 10 Examples of genus- and species-specific probes have been described by Barry et al., 1990, Biotechnology 8:233-236, Atlas and Bej, "PCR protocols: A guide to method and application," p. 399-406; and in Gen-probe international patent application W088/03957 (these references are incorporated herein by reference). The invention
- 15 claimed in this application differs from these inventions in the range of target detected (all bacteria- fungus-yeast) and the application focus (not clinical).

The use of a panel of rRNA probes, including a universal bacterial probe, gram-positive and gram-negative probes and species or group specific probes provides

20 clinically useful information, not a single universal bacterial probe; since different pathologies, drugs and antibiotic therapy is recommended for various bacterial - fungus-yeast infections.

In these precedent patents, universal bacterial probes are using only for positive controls. Universal primers for bacterial- fungus-yeast are using for identification

25 after cloning and sequencing of the amplified product or hybridization on a DNA chip. The used of rRNA targets for sterility controls for detection of alive bacteria and fungus-yeast in sterile or non sterile industrial product has not been described before this invention.

Preferred universal couple of primers for the one step RT PCR detection of bacteria and fungus-yeast comprise a probing nucleobase sequence selected from the group consisting of

- |    |  |                          |                  |
|----|--|--------------------------|------------------|
|    | Seq ID No 1  | TGGAGCATGTGGTTTAATTCGA   | [primer forward] |
| 5  | Seq ID No 2  | TGCGGGACTTAACCCAACA      | [primer reverse] |
|    | Seq ID No 3  | AGAGTTTGATCATGGCTCAGA    | [primer forward] |
|    | Seq ID No 4  | TTACCCACCTACTAGCTAAT     | [primer reverse] |
|    | Seq ID No 5  | GYGGAGCATGTGGYTAAATTCG   | [primer forward] |
|    | Seq ID No 6  | TTGCGCTCGTTRCGGGACTT     | [primer reverse] |
| 10 | Seq ID No 7  | GGGAAACTCACCAGGTCCA      | [primer forward] |
|    | Seq ID No 8  | CGTTATCGCAATTAAGCAGACA   | [primer reverse] |
|    | Seq ID No 9  | GGTAACGGGGAATWAGGGTTC    | [primer forward] |
|    | Seq ID No 10   | TTGGGTAATTTGCGCGCCTG     | [primer reverse] |
| 15 | Preferred universal probes for the detection of bacteria and fungus-yeast comprise a probing nucleobase sequence selected from the group consisting of |                          |                  |
|    | Seq ID No 11   | TGCATGGYTGTCGTCAGCTCGTG  | [probe forward]  |
|    | Seq ID No 12   | GAGTGGCGGACGGGTGAGTAA    | [probe forward]  |
|    | Seq ID No 13   | ACAGGTGGTGCATGGTTGTC     | [probe forward]  |
| 20 | Seq ID No 14   | TCAGCTCGTGTCGTGAGATGTT   | [probe forward]  |
|    | Seq ID No 15   | ACAGGTGCTGCATGGCTGTC     | [probe forward]  |
|    | Seq ID No 16   | TCAGCTCGTGTTGTGAAATGTT   | [probe forward]  |
|    | Seq ID No 17   | AGGATTGACAGATTGAGAGCTCTT | [probe forward]  |
|    | Seq ID No 18   | CGGAGAGGGAGCCTGAGAA      | [probe forward]  |
| 25 | Seq ID No 19   | CGGCTACCACATCCAAGGAA     | [probe forward]  |

The probes and primer sets, methods and kits of this invention are particularly well suited for use in simplex or multiplex one step RT PCR assays wherein all the

bacteria and/or fungus-yeast in a sample can be detected alive and quantitated. The total number of colony forming units (CFU) of bacteria and/or fungus-yeast can be directly determined.

- 5 The following examples are intended to be illustrative of the various methods and compounds of the invention.

#### **EXAMPLE 1**

- 10 A preferred method for analysis of sample by single filtration (filterable liquids). Specificity of extraction from bacteria or fungus-yeast ribonucleotide from the sample up to 1000 mL by centrifiltration following by an incubation with DNase.

- 1 - The liquid sample (up to 1000 mL) is passed through a polycarbonate membrane  
15 (up to 0,45  $\mu\text{m}$ ) or PVDF membrane (up to 0,45  $\mu\text{m}$ ) or PES membrane (up to 0,45  $\mu\text{m}$ ) via centrifugation (swing rotor) at 2000 g or a vacuum pump.

#### **Enzymatic lysis**

- 20 2 - Transfer the filter in a 50 mL sterile tube with up to 1 mL of enzymatic lysis buffer and then incubate at  $35^{\circ}\text{C} \pm 2^{\circ}\text{C}$  for up to one hour.

- 3 - Centrifuge the liquid (lysis buffer) for 5 minutes at 2000 g.

- 25 **Or mecanic lysis**

- 2 - Recovery bacteria and/or fungus-yeast in 600  $\mu\text{L}$  of media. Incubate for up to one hour at  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ .

3 - Add x mg of acid-washed glass beads (100-500  $\mu\text{m}$  in diameter) and x  $\mu\text{L}$  of lysis buffer. Disrupt cells in the Mixer Mills MM300 (Retsch) for up to 90 minutes at maximal speed.

5

4 - The lysat is processed for RNA purification with commercial kits. Our preferred RNA extraction kit is the " MagNaPure LC RNA isolation kit II " on the workstation MagNaPure LC<sup>TM</sup> (Roche Diagnostics). The elution volume is up to 100  $\mu\text{L}$ . Incubation with DNase is processed during the purification.

10

5 – 2  $\mu\text{L}$  (up to 5  $\mu\text{L}$ ) of pure RNA extract is used for the one step real time RT-PCR (LightCycler<sup>TM</sup> ) with enzyme like Tth and the following program with Taqman Probe :

15	I : Reverse transcription	61°C / 20 min (20°C/sec)
	II : Denaturation	95°C / 30 secondes (20°C/sec)
	III : PCR (35 cycles)	95°C / 5 seconds (20°C/sec)
		60°C / 30 seconds (20°C/sec)

The emitted fluorescence is measured at the end of the 60 seconds.

20

6 – 2  $\mu\text{L}$  (up to 5  $\mu\text{L}$ ) of pure RNA extract is used for the one step real time RT-PCR (LightCycler<sup>TM</sup> ) with enzyme like Tth and the following program with Hybridization Probe :

25	I : Reverse transcription	61°C / 20 min (20°C/sec)
	II : Denaturation	95°C / 30 secondes (20°C/sec)
	III : PCR (35 cycles)	95°C / 2 seconds (20°C/sec)
		58°C / 8 seconds (20°C/sec)



16

72°C / 16 secondes (20°C/sec)

The emitted fluorescence is measured at the end of the 8 seconds.

## EXAMPLE 2

5

A preferred method for analysis of sample by centrifugation ( non filterable liquids). Specificity of extraction from bacteria or fungus-yeast ribonucleotide from the sample up to 1000 mL by centrifugation.

10

### Enzymatic lysis

1 – Centrifuge the sample for 5 minutes at 11000 g, discard the supernatant.

15

2 – Resuspend pellet in x  $\mu$ L of lysis buffer. Vortex and incubate for up to one hour at 37°C  $\pm$  1°C.

3 – Add x  $\mu$ L of lysis buffer. Vortex and incubate for up to 30 minutes at 50°C  $\pm$  2°C .

20

### Or mecanic lysis

1 – Centrifuge the sample for 5 minutes at 11000 g, discard the supernatant.

2 – Resuspend pellet in x  $\mu$ L of lysis buffer.

25

3 - Add x mg of acid-washed glass beads (100-500  $\mu$ m in diameter). Disrupt cells in the Mixer Mills MM300 (Retsch) for up to 90 minutes at maximal speed.

4 - The lysat is processed for RNA purification with commercial kits. Our preferred RNA extraction kit is the " MagNaPure LC RNA isolation kit II " on the workstation MagNaPure LC™ (Roche Diagnostics). The elution volume is up to 100 µL.

5 Incubation with DNase is processed during the purification.

5 – 2 µL (up to 5 µL) of pure RNA extract is used for the one step real time RT-PCR (LightCycler™ ) with enzyme like Tth and the following program with Taqman Probe :

10

I : Reverse transcription	61°C / 20 min (20°C/sec)
II : Denaturation	95°C / 30 secondes (20°C/sec)
III : PCR (35 cycles)	95°C / 5 seconds (20°C/sec)
	60°C / 30 seconds (20°C/sec)

15 The emitted fluorescence is measured at the end of the 60 seconds.

6 – 2 µL (up to 5 µL) of pure RNA extract is used for the one step real time RT-PCR (LightCycler™ ) with enzyme like Tth and the following program with Hybridization Probe :

20

I : Reverse transcription	61°C / 20 min (20°C/sec)
II : Denaturation	95°C / 30 secondes (20°C/sec)
III : PCR (35 cycles)	95°C / 2 seconds (20°C/sec)
	58°C / 8 seconds (20°C/sec)
	72°C / 16 secondes (20°C/sec)

25

The emitted fluorescence is measured at the end of the 8 seconds.

**EXAMPLE 3**

A preferred method for analysis of sample by direct lysis and recovery of RNA on DiEthylAminoEthyl cellulose DEAE membrane ( non filterable liquids).

- 5    Specificity of extraction from bacteria or fungus-yeast ribonucleotide from the sample up to 1000 mL by centrifugation on DEAE membrane of the lysat following by an incubation with DNase.

**Enzymatic lysis**

10

1 - Add x mL of lysis buffer in the sample. Vortex and incubate for up to one hour at  $35^{\circ}\text{C} \pm 2^{\circ}\text{C}$ .

15

2 - Add x mL of lysis buffer. Vortex and incubate for up to 30 minutes at  $50^{\circ}\text{C} \pm 2^{\circ}\text{C}$ .

**Or mecanic lysis**

20

1 - Add x mg of acid-washed glass beads (100-500  $\mu\text{m}$  in diameter) and x  $\mu\text{L}$  of lysis buffer.

2 - Disrupt cells in the Mixer Mills MM300 (Retsch) for up to 90 minutes at maximal speed.

- 25    3 – After complete lysis, add the lysat on the pre-filter membrane (polypropylene 10  $\mu\text{m}$ ) to retain particules with a size superior at 10  $\mu\text{m}$ . The liquid filtered is charged in nucleic acid.

4 - Add the the lysat on the pre-wased DEAE membrane to retain all nucleic acid.  
Centrifuge and wash the membrane.

5 - Add salt (up to 1 mL) to recovery nucleic acid in a sterile clean tube by  
5 centrifugation or with a vacuum pump.

6 - The lysat is processed for RNA purification with commercial kits. Our preferred  
RNA extraction kit is the " MagNaPure LC RNA isolation kit II " on the workstation  
MagNaPure LC™ (Roche Diagnostics). The elution volume is up to 100 µL.  
10 Incubation with DNase is processed during the purification.

7 – 2 µL (up to 5 µL) of pure RNA extract is used for the one step real time RT-PCR  
(LightCycler™ ) with enzyme like Tth and the following program with Taqman  
Probe :

15

I : Reverse transcription	61°C / 20 min (20°C/sec)
II : Denaturation	95°C / 30 secondes (20°C/sec)
III : PCR (35 cycles)	95°C / 5 seconds (20°C/sec)
	60°C / 30 seconds (20°C/sec)

20 The emitted fluorescence is measured at the end of the 60 seconds.

8 – 2 µL (up to 5 µL) of pure RNA extract is used for the one step real time RT-PCR  
(LightCycler™ ) with enzyme like Tth and the following program with Hybridization  
Probe :

25

I : Reverse transcription	61°C / 20 min (20°C/sec)
II : Denaturation	95°C / 30 secondes (20°C/sec)
III : PCR (35 cycles)	95°C / 2 seconds (20°C/sec)

20

58°C / 8 seconds (20°C/sec)

72°C / 16 secondes (20°C/sec)

The emitted fluorescence is measured at the end of the 8 seconds.

## CLAIMS

- 5 1. A method and kit for determining the presence of bacteria or fungus-yeast ribonucleic acid (RNA) in a sample suspected of containing said bacteria and/or fungus, wherein said polynucleotide comprises a selected target region, said method comprising:
- 10 (a) extract bacteria or fungus-yeast ribonucleic acid (RNA) from the sample up to 1000 ml by centrifiltration on membranes and /or DEAE resin following by incubation with DNase.
- (b) incubating the bacteria or fungus-yeast ribonucleic acid (RNA) with a  
15 thermostable enzyme with RNA-dependent Reverse Transcriptase activity and with DNA-dependent Polymerase activity, allowing the combination of RT and PCR in a single-tube reaction, such as Tth DNA polymerase, and polynucleotide primers with a nucleotide sequence selected from the group consisting of
- |    |              |                        |                  |
|----|--------------|------------------------|------------------|
| 20 | Seq ID No 2  | TGCGGGACTTAACCCAACA    | [primer reverse] |
|    | Seq ID No 4  | TTACCCACCTACTAGCTAAT   | [primer reverse] |
|    | Seq ID No 6  | TTGCGCTCGTTRCGGGACTT   | [primer reverse] |
|    | Seq ID No 8  | CGTTATCGCAATTAAGCAGACA | [primer reverse] |
|    | Seq ID No 10 | TTGGGTAATTTGCGCGCCTG   | [primer reverse] |
- 25 under conditions which allow hybridization of the polynucleotide to the ribonucleotide target region and Reverse Transcriptase activity of said DNA polymerase for cDNA synthesis; and

(c) amplified the cDNAs formed to a detectable level by Polymerase Chain Reaction with said DNA polymerase and polynucleotide primers and probes with a nucleotide sequence selected from the group consisting of

5	Seq ID No 1	TGGAGCATGTGGTTTAATTCGA	[primer forward]
	Seq ID No 2	TGCGGGACTTAACCCAACA	[primer reverse]
	Seq ID No 3	AGAGTTTGATCATGGCTCAGA	[primer forward]
	Seq ID No 4	TTACCCACCTACTAGCTAAT	[primer reverse]
	Seq ID No 5	GYGGAGCATGTGGYTAAATTCG	[primer forward]
10	Seq ID No 6	TTGCGCTCGTTRCGGGACTT	[primer reverse]
	Seq ID No 7	GGGAACTCACCAGGTCCA	[primer forward]
	Seq ID No 8	CGTTATCGCAATTAAGCAGACA	[primer reverse]
	Seq ID No 9	GGTAACGGGGAATWAGGGTTC	[primer forward]
	Seq ID No 10	TTGGGTAATTTGCGCGCCTG	[primer reverse]
15	Seq ID No 11	TGCATGGYTGTCTCAGCTCGTG	[probe forward]
	Seq ID No 12	GAGTGGCGGACGGGTGAGTAA	[probe forward]
	Seq ID No 13	ACAGGTGGTGCATGGTTGTC	[probe forward]
	Seq ID No 14	TCAGCTCGTGTCGTGAGATGTT	[probe forward]
	Seq ID No 15	ACAGGTGCTGCATGGCTGTC	[probe forward]
20	Seq ID No 16	TCAGCTCGTGTTGTGAAATGTT	[probe forward]
	Seq ID No 17	AGGATTGACAGATTGAGAGCTCTT	[probe forward]
	Seq ID No 18	CGGAGAGGGAGCCTGAGAA	[probe forward]
	Seq ID No 19	CGGCTACCACATCCAAGGAA	[probe forward]

- 25 2. The method and kit of claim 1, wherein the cDNA target sequence synthesised by Reverse Transcriptase activity of the enzyme like Tth polymerase is amplified by the DNA-dependent Polymerase activity of DNA polymerase in the same tube by means of one step real time RT-PCR.

3. The method and kit of claim 1 and 2, wherein the composition for detecting bacteria comprising a polynucleotide primers and a probe consisting of the sequence

Seq ID No 1 TGGAGCATGTGGTTTAATTCGA [primer forward]

5 Seq ID No 2 TGCGGGACTTAACCCAACA [primer reverse]

Seq ID No 11 TGCATGGYTGTCTGTCAGCTCGTG [probe forward]

4. The method and kit of claim 1 and 2, wherein the composition for detecting bacteria comprising a polynucleotide primers and a probe consisting of the sequence

10 Seq ID No 3 AGAGTTTGATCATGGCTCAGA [primer forward]

Seq ID No 4 TTACCCACCTACTAGCTAAT [primer reverse]

Seq ID No 12 GAGTGGCGGACGGGTGAGTAA [probe forward]

5. The method and kit of claim 1 and 2, wherein the composition for detecting bacteria comprising a polynucleotide primers and a probe consisting of the sequence

15 Seq ID No 5 GYGGAGCATGTGGYTAAATTCG [primer forward]

Seq ID No 6 TTGCGCTCGTTRCGGGACTT [primer reverse]

Seq ID No 13 ACAGGTGGTGCATGGTTGTC [probe forward]

Seq ID No 14 TCAGCTCGTGTCGTGAGATGTT [probe forward]

20 Seq ID No 15 ACAGGTGCTGCATGGCTGTC [probe forward]

Seq ID No 16 TCAGCTCGTGTTGTGAAATGTT [probe forward]

6. The method and kit of claim 1 and 2, wherein the composition for detecting fungus-yeast comprising a polynucleotide primers and a probe consisting of the sequence

25 Seq ID No 7 GGGAAACTCACCAGGTCCA [primer forward]

Seq ID No 8 CGTTATCGCAATTAAGCAGACA [primer reverse]

Seq ID No 17 AGGATTGACAGATTGAGAGCTCTT [probe forward]



7. The method and kit of claim 1 and 2, wherein the composition for detecting fungus-yeast comprising a polynucleotide primers and a probe consisting of the sequence

5	Seq ID No 9 GGTAACGGGGAATWAGGGTTC	[primer forward]
	Seq ID No 10 TTGGGTAATTTGCGCGCCTG	[primer reverse]
	Seq ID No 18 CGGAGAGGGAGCCTGAGAA	[probe forward]
	Seq ID No 19 CGGCTACCACATCCAAGGAA	[probe forward]

10 8. The method and kit of one of claims 1 to 6, wherein the preferred combination of primers and probes used for detection all bacteria and/or fungus-yeast consisting of the sequence :

Seq ID No 1+ Seq ID No 2 +Seq ID No 11

or

15 Seq ID No 3+ Seq ID No 4 +Seq ID No 12

or

Seq ID No 5+ Seq ID No 6 +Seq ID No 13 + Seq ID No 14 + Seq ID No 15 +Seq ID No 16

or

20 Seq ID No 7+ Seq ID No 8 +Seq ID No 17

or

Seq ID No 9+ Seq ID No 10 +Seq ID No 18 + Seq ID No 19

or

Seq ID No 1+ Seq ID No 2 +Seq ID No 11 + Seq ID No 7+ Seq ID No 8 +Seq ID No

25 17

or

Seq ID No 3+ Seq ID No 4 +Seq ID No 12 + Seq ID No 7+ Seq ID No 8 +Seq ID No

17

or

Seq ID No 5+ Seq ID No 6 +Seq ID No 13 + Seq ID No 14 + Seq ID No 15 +Seq ID  
No 16 + Seq ID No 9+ Seq ID No 10 +Seq ID No 18 + Seq ID No 19

- 5     9. The method and kit of one of claims 1 to 8, wherein the polynucleotide primers  
and probes are natural nucleic acid or Peptide Nucleic Acid (PNA) which can  
hybridize to nucleic acid (DNA and RNA).
- 10    10. The method and kit of one of claims 1 to 9, and also quantified this RNA for a  
comparison with quantified external standard RNA from by exemple *Escherichia*  
*coli* and *Candida spp.*

## SEQUENCE LISTING

<110> Genolife

<120> One step real-time RT PCR kits for the universal detection of organisms in industrial products

<130> D20782

<150> US 60/425,327

<151> 2002-11-12

<160> 19

<170> PatentIn version 3.2

<210> 1

<211> 22

<212> DNA

<213> Artificial

<220>

<221> Source

<222> (1)..(22)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 1

tggagcatgt ggtttaattc ga

22

<210> 2

<211> 19

<212> DNA

<213> Artificial

<220>

<221> source

<222> (1)..(19)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 2

tgcgggactt aaccaaca

19

<210> 3

<211> 21

<212> DNA

<213> Artificial

<220>

<221> source

<222> (1)..(21)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 3

agagtttgat catggctcag a

21

<210> 4  
<211> 21  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(21)  
<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 4  
ttaccccacc tactagctaa t 21

<210> 5  
<211> 22  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(22)  
<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 5  
gyggagcatg tggyttaatt cg 22

<210> 6  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(20)  
<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 6  
ttgcgctcgt trcgggactt 20

<210> 7  
<211> 19  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(19)  
<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 7  
gggaaactca ccagggtcca 19

<210> 8

<211> 22  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(22)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 8  
cgttatcgca attaagcaga ca 22

<210> 9  
<211> 21  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(21)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 9  
ggtaacgggg aatwagggtt c 21

<210> 10  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(20)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 10  
ttgggtaatt tgcgcgcctg 20

<210> 11  
<211> 23  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(23)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 11  
tgcattggytg tcgtcagctc gtg 23

<210> 12  
<211> 21

<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(21)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 12  
gagtggcgga cgggtgagta a 21

<210> 13  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(20)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 13  
acaggtggtg catggttgtc 20

<210> 14  
<211> 22  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(22)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 14  
tcagctcgtg tcgtgagatg tt 22

<210> 15  
<211> 20  
<212> DNA  
<213> Artificial

<220>  
<221> source  
<222> (1)..(20)  
<223> /note = "Description of artificial sequence: oligonucleotide  
primer"

<400> 15  
acaggtgctg catggctgctc 20

<210> 16  
<211> 22  
<212> DNA

<213> Artificial

<220>

<221> source

<222> (1) .. (22)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 16

tcagctcgtg ttgtgaaatg tt

22

<210> 17

<211> 24

<212> DNA

<213> Artificial

<220>

<221> source

<222> (1) .. (24)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 17

aggattgaca gattgagagc tctt

24

<210> 18

<211> 19

<212> DNA

<213> Artificial

<220>

<221> source

<222> (1) .. (19)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 18

cgagagaggga gcctgagaa

19

<210> 19

<211> 20

<212> DNA

<213> Artificial

<220>

<221> source

<222> (1) .. (20)

<223> /note = "Description of artificial sequence: oligonucleotide primer"

<400> 19

cggtaccac atccaaggaa

20